BD BACTEC™ Plus Aerobic/F Culture Vials and BACTEC™ Plus Anaerobic/F Culture Vials

Soybean-Casein Digest Broth

English: pages 1 – 4 Español: páginas 13 – 16
 Français : pages 4 – 7 Dansk: side 16 – 19
 Deutsch: Seiten 7 – 10 Português: páginas 19 – 22
 Italiano: pagine 10 – 13 Svenska: sidan 22 – 25

8085859(02)
 2015-04

Contact your local BD representative for instructions. / Свържете се с местния представител на BD за инструкции. / Pokyny vám poskytne místní zástupce společnosti BD. / Kontakt den lokale BD repræsentant for at få instruktioner. / Επικοινωνήστε με τον τοπικό αντιπρόσωπο της BD για οδηγίες. / Kasutusjuhiste suhtes kontakteeruge oma kohaliku BD esindajaga. / Ota yhteyks lähimpään BD:n edustajaan ohjeiden saamiseksi. / Kontaktiraj lokalnog predstavnika BD za upute. / A használati utasítást kérje a BD helyi képviselőjétől. / Нұсқаулар үшін жергілікті BD өкілімен хабарласыңыз. / Lai saņemtu norādījumus, sazinieties ar vietējo BD pārstāvi. / Naudojimo instrukcijų teiraukitės vietos BD įgaliotojo atstovo. / Neem contact op met uw plaatselijke BD-vertegenwoordiger voor instructies. / Kontakt din lokale BD-representant for mer informasjon. / Aby uzyskać instrukcje użytkowania, skontaktuj się z lokalnym przedstawicielem BD. / Contacte o representante local da BD para instruções. / Pentru instrucțiuni, contactați reprezentantul local BD. / Для получения указаний обратитесь к местному представителю компании BD. / Instrukcie získate u miestneho zástupcu spoločnosti BD. / Obratite se svom lokalnom predstavniku kompanije BD za uputstva. / Kontakta närmaste BD-representant för anvisningar. / Talimatlar için yerel BD temsilcinize temasa geçin. / За інструкціями зверніться до місцевого представника компанії BD.

INTENDED USE

BD BACTEC Plus Aerobic/F and Plus Anaerobic/F media are used in a qualitative procedure for the aerobic and anaerobic culture and recovery of microorganisms (bacteria and yeast) from blood. The principal use of these media is with the BD BACTEC fluorescent series instruments.

SUMMARY AND EXPLANATION

The sample to be tested is inoculated into one or more vials which are inserted into the BACTEC fluorescent series instrument for incubation and periodic reading. Each vial contains a chemical sensor which can detect increases in CO₂ produced by the growth of microorganisms. The sensor is monitored by the instrument every ten minutes for an increase in its fluorescence, which is proportional to the amount of CO₂ present. A positive reading indicates the presumptive presence of viable microorganisms in the vial. Detection is limited to microorganisms that will grow in a particular type of medium.

Resins have been described for the treatment of blood specimens both prior to and after their inoculation into culture media. Resins have been incorporated into BACTEC culture media to enhance recovery of organisms without a need for special processing.¹⁻³

PRINCIPLES OF THE PROCEDURE

If microorganisms are present in the test sample inoculated into the BACTEC vial, CO₂ will be produced when the organisms metabolize the substrates present in the vial. Increases in the fluorescence of the vial sensor caused by the higher amount of CO₂ are monitored by the BACTEC fluorescent series instrument. Analysis of the rate and amount of CO₂ increase enables the BACTEC fluorescent series instrument to determine if the vial is positive, i.e., that the test sample contains viable organisms.

REAGENTS

The BACTEC culture vials contain the following reactive ingredients prior to processing:

List of Ingredients	BACTEC Plus Aerobic/F	BACTEC Plus Anaerobic/F
Processed Water	30 mL* w/v	25 mL w/v
Soybean-Casein Digest Broth	3.0%	3.0%
Yeast Extract	0.25%	0.4%
Animal Tissue Digest	—	0.05%
Amino Acids	0.05%	0.25%
Sugar	—	0.25%
Sodium Citrate	—	0.02%
Sodium Polyanetholsulfonate (SPS)	0.05%	0.05%
Vitamins	0.025%	0.0006%
Antioxidants/Reductants	0.005%	0.16%
Nonionic Adsorbing Resin	16.0%	16.0%
Cationic Exchange Resin	1.0%	1.0%

All BACTEC media are dispensed with added CO₂. Anaerobic media are prereduced and dispensed with CO₂ and N₂.

*The medium volume has been increased in the BACTEC Plus Aerobic/F from 25 mL to 30 mL.

Warnings and Precautions:

The prepared culture vials are for *in vitro* diagnostic use.

This Product Contains Dry Natural Rubber.

Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"⁴⁻⁷ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.

Prior to use, each vial should be examined for evidence of contamination such as cloudiness, bulging or depressed septum, or leakage. **DO NOT USE** any vial showing evidence of contamination. A contaminated vial could contain positive pressure. If a contaminated vial is used for direct draw, gas or contaminated culture media could be refluxed into the patient's vein. Vial contamination may not be readily apparent. If a direct draw procedure is used, monitor the process closely to avoid refluxing materials into the patient.

Prior to use, the user should examine the vials for evidence of damage or deterioration. Vials displaying turbidity, contamination, or discoloration (darkening) should not be used. On rare occasions, the glass bottle neck may be cracked and the neck may break during removal of the flip-off cap or in handling. Also, on rare occasions a vial may not be sealed sufficiently. In both cases the contents of the vials may leak or spill, especially if the vial is inverted. If the vial has been inoculated, treat the leak or spill with caution, as pathogenic organisms/agents may be present. Before discarding, sterilize all inoculated vials by autoclaving.

Positive culture vials for subculturing or staining, etc.: before sampling it is necessary to release gas which often builds up due to microbial metabolism. Sampling should be performed in a biological safety cabinet if possible, and appropriate protective clothing, including gloves and masks, should be worn. See Procedure section for more information on subculturing.

To minimize the potential of leakage during inoculation of specimen into culture vials, use syringes with permanently attached needles or **Luer-Lok™** tips.

Storage Instructions

The **BACTEC** vials are ready for use as received and require no reconstitution or dilution. Store in a cool, dry place (2° – 25 °C), **out of direct light**.

SPECIMEN COLLECTION

The specimen must be collected using sterile techniques to reduce the chance of contamination. The recommended specimen volume is 8 – 10 mL. It is recommended that the specimen be inoculated into the **BACTEC** vials at bedside. A 10cc or 20cc syringe with a **Luer-Lok** brand tip is used to draw the sample, or a **Vacutainer™** Brand Needle Holder and a **Vacutainer™** Brand Blood Collection Set, **Vacutainer™ Safety-Lok™** Blood Collection Set or other tubing "butterfly" set may be used. If using a needle and tubing set (direct draw), carefully observe the direction of blood flow when starting sample collection. The vacuum in the vial will usually exceed 10 mL, so the user should monitor the volume collected by means of the 5 mL graduation marks on the vial label. Sample volumes as low as 3 mL can be used, however, recovery will not be as great as with larger volumes. **The inoculated BACTEC vial should be transported to the laboratory as quickly as possible.**

PROCEDURE

Remove the flip-off cap from the **BACTEC** vial top and inspect the vial for cracks, contamination, excessive cloudiness, and bulging or indented septums. **DO NOT USE** if any defect is noted. Before inoculating, swab the septum with alcohol (iodine is **not** recommended). Aseptically inject or draw directly 8 – 10 mL of specimen per vial. If sample volumes of 3 – 7 mL are used, recovery will not be as great as with larger volumes (see Limitations of the Procedure). **Inoculated aerobic and anaerobic vials should be placed in the BACTEC fluorescent series instrument as soon as possible** for incubation and monitoring. If placement of an inoculated vial into the instrument has been delayed and visible growth is apparent, it should not be tested in the **BACTEC** fluorescent series instrument, but rather it should be subcultured, Gram-stained and treated as a presumptively positive bottle.

Vials entered into the instrument will be automatically tested every ten minutes for the duration of the testing protocol period. Positive vials will be determined by the **BACTEC** fluorescent series instrument and identified as such (see the appropriate **BACTEC** Fluorescent Series Instrument User's Manual). The sensor inside the bottle will not appear visibly different in positive and negative vials, however the **BACTEC** fluorescent series instrument can determine a difference in fluorescence.

If at the end of the testing period a negative vial appears visually positive (i.e., chocolateized blood, bulging septum, lysed and/or very darkened blood in **BACTEC** Plus Aerobic medium), it should be subcultured and Gram-stained and treated as a presumptive positive.

Positive vials should be subcultured and Gram-stained. In a great majority of cases, organisms will be seen and a preliminary report can be made to the physician. Subcultures to solid media and a preliminary direct antimicrobial susceptibility test may be prepared from fluid in the **BACTEC** vials.

Subculturing: Prior to subculturing, put the vial in an upright position, and place an alcohol wipe over the septum. To release pressure in the vial, use an appropriate venting unit (BD catalog # 249560) or equivalent. The needle should be removed after the pressure is released and before sampling the vial for subculture. The insertion and withdrawal of the needle should be done in a straight-line motion, avoiding any twisting motions.

QUALITY CONTROL

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent CLSI guidance and CLIA regulations for appropriate Quality Control practices.

DO NOT USE culture vials past their expiration date.

DO NOT USE culture vials that exhibit any cracks or defects; discard the vial in the appropriate manner.

Quality Control Certificates are provided with each carton of media. Quality Control Certificates list test organisms, including ATCC® cultures specified in the CLSI Standard M22, *Quality Control for Commercially Prepared Microbiological Culture Media*. The range of time-to-detection in hours was ≤ 72 hours for each of the organisms listed on the Quality Control Certificate for this medium:

Aerobic Medium Organisms

<i>Neisseria meningitidis</i> ATCC 13090	<i>Candida glabrata</i> ATCC 66032
<i>Haemophilus influenzae</i> ATCC 19418	<i>Staphylococcus aureus</i> ATCC 25923
<i>Streptococcus pneumoniae</i> * ATCC 6305	<i>Escherichia coli</i> ATCC 25922
<i>Streptococcus pyogenes</i> ATCC 19615	<i>Alcaligenes faecalis</i> ATCC 8750
<i>Pseudomonas aeruginosa</i> ATCC 27853	

Anaerobic Medium Organisms

<i>Clostridium perfringens</i> ATCC 13124	<i>Streptococcus pneumoniae</i> ATCC 6305
<i>Bacteroides fragilis</i> * ATCC 25285	<i>Escherichia coli</i> ATCC 25922
<i>Bacteroides vulgatus</i> ATCC 8482	<i>Staphylococcus aureus</i> ATCC 25923

*CLSI recommended strain

For information on Quality Control for the **BACTEC** fluorescent series instrument, refer to the appropriate **BACTEC** Fluorescent Series Instrument User's Manual.

RESULTS

A positive sample is determined by the **BACTEC** fluorescent series instrument and indicates the presumptive presence of viable microorganisms in the vial.

LIMITATIONS OF THE PROCEDURE

Contamination

Care must be taken to prevent contamination of the sample during collection and inoculation into the **BACTEC** vial. A contaminated sample will give a positive reading, but will not indicate a relevant clinical sample. Such a determination must be made by the user based on such factors as type of organism recovered, occurrence of the same organism in multiple cultures, patient history, etc.

Recovery of SPS Sensitive Organisms From Blood Samples

Because blood can neutralize the toxicity of SPS toward organisms sensitive to SPS, the presence of maximum volumes of blood (8 – 10 mL) can help to optimize recovery of these organisms. To enhance the growth of SPS sensitive organisms when less than 8 mL of blood is inoculated, additional whole human blood may be added.

Some fastidious organisms, such as certain *Haemophilus* species, require growth factors, such as NAD, or factor V, which are provided by the blood specimen. If the blood specimen volume is 3.0 mL or less, an appropriate supplement may be required for recovery of these organisms. **BACTEC FOS™** Fastidious Organism Supplement or whole human blood may be used as nutritional supplements.

Nonviable Organisms

A Gram-stained smear from a culture medium may contain small numbers of nonviable organisms derived from media constituents, staining reagents, immersion oil, glass slides, and specimens used for inoculation. In addition, the patient specimen may contain organisms that will not grow in the culture medium or in media used for subculture. Such specimens should be subcultured to special media as appropriate.

Antimicrobial Activity

Neutralization of the antimicrobial activity by resins varies depending upon dosage level and timing of specimen collection. The use of supplementary additives should be considered in appropriate situations; as an example, the addition of penicillinase when β -lactam therapy is being employed.

For information on antimicrobial agents neutralized by **BACTEC** resins, contact the BD Technical Services department at the toll free number listed below.

Recovery of *Streptococcus pneumoniae*

In aerobic media, *S. pneumoniae* will typically be visually and instrument positive, but in some cases no organism will be seen on Gram stain or recovered on routine subculture. If an anaerobic vial was also inoculated, the organism can usually be recovered by performing an aerobic subculture of the anaerobic vial, since this organism has been reported to grow well under anaerobic conditions.¹¹

General Considerations

Optimum recovery of isolates will be achieved by adding maximum amounts of blood. Use of lower volumes may adversely affect recovery and/or detection times of organisms such as *Peptostreptococcus*, *Peptococcus*, *Bacteroides asaccharolyticus*. Blood may contain antimicrobials or other inhibitors which may slow or prevent the growth of microorganisms. False negative readings may result when certain organisms are present which do not produce enough CO₂ to be detected by the system or if significant growth has occurred before placing the vial into the system. False positivity may occur when the white blood cell count is high.

Due to the nature of biological materials in media products and inherent organism variability, the user should be cognizant of potential variable results in the recovery of certain microorganisms.

EXPECTED VALUES AND SPECIFIC PERFORMANCE CHARACTERISTICS

Performance of the **BD BACTEC Plus Aerobic/F** and **Plus Anaerobic/F** media has been established by a number of external clinical studies.^{1-3,8,9} Seeded laboratory studies performed by BD have shown equivalent performance of the **BD BACTEC Plus Aerobic/F** and **Plus Anaerobic/F** media to nonradiometric **BACTEC Plus 26** and **Plus 27**.¹⁰

AVAILABILITY

Cat. No.	Description
442192	BD BACTEC™ Plus Aerobic/F Medium
442193	BD BACTEC™ Plus Anaerobic/F Medium

REFERENCES

1. Wallis, C. et al. 1980. Rapid isolation of bacteria from septicemic patients by use of an antimicrobial agent removal device. *J. Clin. Microbiol.* 11:462-464.
 2. Applebaum, P.C. et al. 1983. Enhanced detection of bacteremia with a new **BACTEC** resin blood culture medium. *J. Clin. Microbiol.* 17:48-51.
 3. Pohlman, J.K. et al. 1995. Controlled clinical comparison of Isolator and **BACTEC** 9240 Aerobic/F resin bottle for detection of bloodstream infections. *J. Clin. Microbiol.* 33:2525-2529.
 4. Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
 5. Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.
 6. U.S. Department of Health and Human Services. 2007. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 5th ed. U.S. Government Printing Office, Washington, D.C.
 7. Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
 8. Smith, J.A. et al. 1995. Comparison of **BACTEC** 9240 and BacT/Alert blood culture systems in an adult hospital. *J. Clin. Microbiol.* 33:1905-1908.
 9. Flayhart, D. et al. 2007. Comparison of **BACTEC** Plus blood culture media to BacT/Alert FA blood culture media for detection of bacterial pathogens in samples containing therapeutic levels of antibiotics. *J. Clin. Microbiol.* 45:816-821.
 10. Data available from BD Diagnostics.
 11. Howden, R.J. 1976. Use of anaerobic culture for the improved isolation of *Streptococcus pneumoniae*. *J. Clin. Pathol.* 29:50-53.
- Technical Information: In the United States contact BD Technical Service and Support at 800-638-8663 or www.bd.com/ds.

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Becton, Dickinson and Company
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113
Australia

ATCC is a trademark of the American Type Culture Collection.

BD, BD Logo, and all other trademarks are property of Becton, Dickinson and Company. © 2015 BD.